

**Tomorrow's leading
technology.
Available today.**



SilaMill Z4

Precise. Fast. Independent. Economical.



Highest Precision

- Milling and grinding in ultra HD
- Heavy industrial quality
- 3 microns repetition accuracy



Fastest Production

- Restorations in under 10 minutes
- 2-second block insertion
- 100,000 RPM electrical high frequency spindle



Complete Independence

- 30 block materials from 12 manufacturers – and growing
- >140 titanium and CoCr prefabricated abutments from 11 manufacturers
- Validated with all established scanners and design software
- Integrated PC with touch screen and Wi-Fi – no laptop or tablet necessary
- Built-in compressed air – no compressor needed



Extremely Economical

- Clear water – no additives necessary
- Easy to learn, easy to operate
- The perfect complement, made by SILADENT: long-life precision tools with competitive pricing



Made in Germany

 **SILADENT**